

The Rexton Day family of hearing instruments features enough technology and flexibility to suit a wide range of hearing needs. Perfect for both new and experienced wearers, Day combines the right technology and performance with an affordable price so you can be confident you are getting the best value possible.

#### ITE AND BTE FEATURES

- 6 Channels
- 3 Programs
- Wireless
- Noise Reduction
- Adaptive Directional Microphone
- Adaptive Feedback Cancellation
- AGC-O control
- DataLogging
- Power-on delay
- Low battery beeps
- Program change beeps

#### BTE ONLY FEATURES:

- Push Button  
(DAY FREE, P AND HP ONLY)
- Volume Control  
(DAY P/HP ONLY)
- Telecoil (DAY FREE, P, HP)
- Battery On/Off Switch
- Battery Door Lock  
(DAY P/HP ONLY)
- Available Colors: Beige, granite, grey, brown, silver

#### CUSTOM OPTIONS

- Push button
- Wireless
- Maximum 4 programmable memories (push button)
- Programmable volume control
- Automatic and adaptive directional microphone
- Programmable telecoil

#### ACCESSORIES

- Remote control unit



*Day*  
Nourish your desire to hear better.



*Day*  
Welcome to a New Day in Better Hearing.



Rev. 3/09 PN 104250374/a

**REXTON**

**REXTON**

Unlock a new Day of better hearing.



## A New Day in Better Hearing

Better hearing doesn't have to be complicated or expensive. The day you accept your hearing loss is the day you can begin to enjoy life again with an appropriate hearing system. If you are looking for an easy-to-use, comfort-focused solution, Rexton Day hearing instruments could be the answer for you.

## Oh Happy Day

Rexton strives to help improve the lives of people with hearing loss using advanced digital processing, fitting flexibility and reliability. Day is a perfect example of our commitment to understand today's value-conscious wearer. Better hearing, better value, and greater satisfaction are the keys to Rexton's 50 years of success. Get started today.

## Hear Better Every Day

Today is the perfect day to realize better hearing doesn't have to be complicated or expensive. Day hearing instruments come in many different styles and are the perfect choice for many individuals. Make your hearing experience the best it can be with any easy to wear Rexton Day hearing solution.



All Day digital circuits are available in a full range of cosmetically appealing models. From the Completely-in-the-Canal (CIC) to Open Fit and Behind-the-Ear (BTE) instruments. Ask your hearing professional which one is the right fit for you.



OPEN BTE

FREE BTE

P BTE

HP BTE

## Fill up your Day with features

### Automatic Adaptive Directional Microphone System

Lively gatherings with friends and conversations in noise are often very difficult for hearing instrument wearers. Thanks to Day's advanced directional microphone system, you can enjoy easier listening in any situation.

### Noise Reduction

Whether you're working in a busy office, waiting for a train in a noisy city subway station, or enjoying a cup of coffee at a sidewalk café, Day's adaptive noise reduction system analyzes incoming signals in real time and isolates important speech signals. Speech is enhanced while distracting background noise is softened, so you can enjoy more effortless listening.

### Adaptive Feedback Cancellation

Day instruments feature an essential Adaptive Feedback Cancellation system. This innovation quickly detects and reduces annoying and embarrassing feedback whistling, before it's heard by you or others nearby.

### DataLogging

By providing a running record of valuable data about your daily listening activities and instrument use, DataLogging helps your hearing care professional provide faster and more accurate fittings. Less time spent adjusting your instrument and fewer follow up visits let you get back to enjoying your day.

### NanoCoating

NanoCoating is a unique technology that helps repel moisture, dirt, and other debris from BTE instruments before it can penetrate their housing. This feature ensures greater durability, fewer repairs, and greater longevity of the instrument



CIC

CANAL

HALF SHELL

FULL SHELL